


# Calcium Kit

## CliniQuant - FSR

### Arsenazo III Method, End Point

 For *in vitro* diagnostic use  
Read this pack insert thoroughly before use

REF	Pack Size	R1 Calcium Arsenazo Reagent	R2 Calcium Standard
CAAFSR-01	2 x 25ml	2 x 25ml	1 x 5ml

#### INTENDED USE

This reagent is intended for quantitative determination of calcium concentration in human serum or plasma and urine.

#### CLINICAL SIGNIFICANCE

In human body, circulating calcium is used for several functions, in skeletal metabolism, neuromuscular function and in hemostasis. Serum calcium is affected by mal-absorption, parathyroidal problems, bone disease, circulating albumin diseases.

#### Typical variations:


Greater than normal levels are seen in hyperparathyroidism, metastatic bone tumor, milk - alkali syndrome, multiple myeloma, Paget's disease, sarcoidosis, vitD intoxication;

Lower than normal levels are seen in hypoparathyroidism, mal-absorption, osteomalacia, pancreatitis, renal failure, vit-D deficiency.

The amount of calcium excreted into the urine reflects intestinal absorption, skeletal resorption, and renal tubular filtration and re-absorption. Under fasting conditions, the intestinal and renal components are relatively fixed and calcium excretion (mg/dl of glomerular filtrate) in the fasting state is used to assess the skeletal components.

#### PRINCIPLE OF THE METHOD

Arsenazo III combines with calcium ions at pH 6.7 to form a blue color complex, the absorbance of which is measured at 620 nm.



#### KIT COMPONENTS

##### Composition

R1 - Calcium Arsenazo Reagent : Arsenazo III 0.167 mmol/l, Imidazole Buffer pH 6.5 99.88 mmol/l

R2 - Calcium Standard : 10 mg/dl, Calcium Carbonate : 0.25 g/l

#### MATERIALS REQUIRED BUT NOT PROVIDED

Laboratory instrumentation, Spectrophotometer UV/VIS with thermostatic cuvette holder or clinical chemistry analyzer: semi automated, calibrated micropipettes, glass or high quality polystyrene cuvettes, test tube/rack, heating bath, controls, saline.

#### REAGENT PREPARATION, STORAGE & STABILITY

Reagent is ready to use. Keep away from direct light sources.

Stability: up to expiration date on labels at 2-8 °C.

#### REAGENT DETERIORATION

1. Discard the reagent if absorbance exceeds 1.2 against distilled water.

2. Keep the Standard vial plugged after use, in order to avoid deterioration.

#### WARNINGS & PRECAUTIONS

1. Reagent may contain some non-reactive and preservative components. It is recommended to handle carefully, avoiding contact with skin and ingestion.

2. Specimens should be considered infectious and handled appropriately.

3. Care must be taken to avoid calcium contamination. The use of plastic tubes or cuvettes free from calcium contamination is strongly recommended. Glassware used be soaked in dilute HCl or dilute HNO<sub>3</sub> and thoroughly rinsed with distilled water.

4. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

#### SPECIMEN

Serum (preferred) plasma (heparinate). Venostasis should be avoided in specimen collection. Samples are stable for 7 days at 2-8°C & 6 months when frozen.

Urine specimen should be collected in 20 to 30 ml of HCl 6M per 24/h specimen (1-2 ml for random urine) in order to prevent calcium salt precipitation.

Dilute urine sample 1:2 with distilled water and multiply results by two.

#### Programme Parameter for MERILYZER CliniQuant

Reading Mode	End Point
Standard Conc.	10 (mg/dl)
Filter - 1 (nm)	620 nm
Temperature	37 °C
Volume (µl)	500
Delay Time (Sec)	5
Reaction Direction	Increase
Reference Low	8.4
Reference High	10.4
Linearity Limit	16

## TEST PROCEDURE

Dispense	Blank	Standard	Sample
Reagent 1	1ml	1ml	1ml
Distilled water	20 µl	-	-
Standard	-	20 µl	-
Sample	-	-	20 µl

Mix, incubate for 1 min at 37°C. Read absorbance of standard (As) and samples (Ax) against reagent blank

## RESULT CALCULATION

Serum/plasma:

Calcium mg/dl = Ax/As x Concentration of Standard

Urine sample:

Calcium mg/dl = Ax/As x Concentration of Standard x 2

24 hours urine sample:

Calcium mg/24h = Ax/As x Concentration of Standard x 2 x urine volume

SI conversion factor: 1 mg/dl x 0.25 = 1 mmol/l

## EXPECTED VALUES

Serum/plasma: 8.4 – 10.4 mg/dl OR 2.1 – 2.6 mmol/l

Urine (men): up to 300 mg/24h OR 7.499 mmol/24h

Urine (women): up to 250 mg/24h OR 6.24 mmol/24h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## QUALITY CONTROL AND CALIBRATION

It is recommended to perform internal quality control with assayed normal (BioNorm) and assayed abnormal (BioPath), to confirm the validity of the test and assure the accuracy of patient result.

Using the recommended Calibrator (BioCal) or the Standard included, calibrate the assay:

- When using a new reagent or lot
- When QC values are out of range

## PERFORMANCE CHARACTERISTICS

### 1. Linearity

The linearity is up to 16 mg/dl or 4 mmol/l.

## 2. Sensitivity/ Limit of detection (LOD)

The limit of detection is 0.04 mg/dl.

The limit of quantification is 0.1 mg/dl.

## 3. Interferences

Gross hemolysis and lipaemia may cause falsely elevated results, a sample blank be set by adding 10 µl sample in 1ml saline.

## 4. Precision

Intra-assay precision

	Mean	SD	CV
n = 20	mg/dl	mg/dl	%
sample 1	7.65	0.02	0.32
sample 2	11.88	0.03	0.29

Inter-assay precision

	Mean	SD	CV
n = 20	mg/dl	mg/dl	%
sample 1	7.89	0.38	4.84
sample 2	12.36	0.58	4.73

## 5. Methods Comparison

Comparison was done between reference Calcium Arsenazo Reagent and CliniQuant - FSR Calcium (A) Reagent (test).

N = 24                      y = 1.204x – 1.021

r<sup>2</sup> = 0.984

## LIMITATIONS

Samples with values above 16 .0mg/dl should be diluted with 0.9% saline, re-run and results multiplied by dilution factor.

## WASTE DISPOSAL

This product is made to be used in professional laboratories. Please consult local regulations for correct waste disposal.

## REFERENCES

- Burtis, C.A., Ashwood, E.R., editors. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia, W.B. Saunders Company, 1994, p. 1887 - 1906.
- Data on file: Meril Diagnostics.

IFU/CAAFSR01/00

06-11-2018

Symbols used on Meril Diagnostics labels:

 REF	Catalogue No.		Attention See Instruction for Use
 LOT	Batch No.		In vitro Diagnostics
	Expiry Date		Consult Instruction for Use
	Manufacturer		Storage Temperature
	Keep Dry		Keep Away from Sunlight
	Manufacturing Date		Do not use if package is damaged
	Authorized European Representative in the European Community		