


Creatine Kinase Kit

CliniQuant - FSR

NAC Activated Method

 For *in vitro* diagnostic use
Read this pack insert thoroughly before use

| REF | Pack Size | R1 CK-NAC Reagent | R2 CK-NAC Reagent |
|-----------|-----------------|-------------------------|-------------------------|
| CKNFSR-01 | 2 x 8 / 2 x 2ml | 2 x 8ml | 2 x 2ml |


INTENDED USE

This reagent is intended for quantitative determination of creatine kinase (CK) level in human serum.

CLINICAL SIGNIFICANCE

Creatine kinase (CK) is an enzyme which is contained in heart, brain and skeletal muscles. Thus, an increase of circulating level of CK may be associated to myocardial infarction, acute cerebrovascular disease, trauma or diseases of skeletal muscles. After a myocardial infarction, CK level begin rising between 4th and 6th hour after first acute symptoms, reaching the peak between 18th and 30th hour and coming back to normal values during the 3rd day.

PRINCIPLE OF THE METHOD



KIT COMPONENTS

Composition :

R1 - CK-NAC Reagent : Imidazole buffer 100mmol/l pH 6.7, AMP 5 mmol/l, A2p5 10 µmol/l, Magnesium acetate 10 mmol/l, ADP 50 mmol/l, NAD 1 mmol/l, NAC 20 mmol/l, Hexokinase > 16 KU/l, G6PDH > 11 Ku/l, Bicine 20mmol/l
R2 - CK-NAC Reagent : Creatine phosphate 263 mol/l, Bicine buffer 20 mmol/l pH 9.3

MATERIALS REQUIRED BUT NOT PROVIDED

Laboratory instrumentation, Spectrophotometer UV/VIS with thermostatic cuvette holder or clinical chemistry analyzer: semi automated, calibrated micropipettes, glass or high quality polystyrene cuvettes, test tube/ rack, heating bath, controls, saline.

REAGENT PREPARATION, STORAGE & STABILITY

Mix reagent 1 & reagent 2 in ratio 4:1. Keep away from direct light sources.

Stability: up to expiration date on labels at 2-8 °C.

Stability of working reagent: 10 days at 2-8 °C.

REAGENT DETERIORATION

1. Discard the reagent if absorbance exceeds 0.6 at 340 nm against distilled water.
2. Discard the working reagent if it fails to achieve assigned assay values of fresh control sera.

WARNINGS AND PRECAUTIONS

1. Do not ingest. Reagent may contain some non-reactive and preservative components. It is recommended to handle carefully, avoiding contact with skin.
2. Specimens should be considered infectious and handled appropriately.
3. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Use unhemolysed serum. Avoid use of plasma containing heparin, EDTA, citrate or flouride. CK activity in serum is unstable and is rapidly lost during storage. Samples are stable for 1 day at 2-8 °C and for longer when frozen at -20°C.

Programme Parameter for MERILYZER CliniQuant

| | |
|--------------------|----------|
| Reading Mode | Rate |
| Factor | 8095 |
| Filter - 1 (nm) | 340 |
| Temperature | 37 °C |
| Volume (µl) | 500 |
| Delay Time (Sec) | 120 |
| Test Time (Sec) | 180 |
| Unit | U/l |
| Reaction Direction | Increase |
| Reference Low | 25 |
| Reference High | 200 |
| Linearity Limit | 1800 |

TEST PROCEDURE

| | |
|---|---------|
| Dispense in tube : Working Reagent | 1000 µl |
| Sample | 20 µl |
| Mix, execute a first reading of absorbance after 2 minutes, incubating at 37°C. Perform other 3 readings at 60 seconds intervals. Calculate the ΔA/min. | |

RESULT CALCULATION

Perform calculations in units per litre, multiplying the ΔA/min by the factor.

Activity in U/l = ΔA/min x 8095

SI conversion factor: 1 U/l x 0.017 = 1 µkat/l

EXPECTED VALUES

Men: 25 – 200 U/l OR 0.4 – 3.5 µkat/l
Women: 25 – 173 U/l OR 0.4 – 3.0 µkat/l

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

QUALITY CONTROL AND CALIBRATION

It is recommended to perform internal quality control with assayed normal (BioNorm) and assayed abnormal (BioPath), to confirm the validity of the test and assure the accuracy of patient result.

When using the recommended Calibrator (BioCal), calibrate the assay:

- When using a new reagent or lot
- When QC values are out of range

PERFORMANCE CHARACTERISTICS

1. Linearity

The linearity is up to 1800 U/l or 31 µkat/l.

2. Sensitivity/ Limit of detection (LOD)

The limit of detection is 5.5 U/l.

The limit of quantification is 16.8 U/l.

3. Interferences

No interference has been observed for the following
Hemoglobin up to 40 mg/dl; Bilirubin up to 40 mg/dl
Triglycerides up to 1000 mg/dl

4. Precision

Intra-assay precision

| | Mean | SD | CV |
|----------|------|------|-----|
| n = 20 | U/l | U/l | % |
| sample 1 | 154 | 0.93 | 0.6 |
| sample 2 | 488 | 1.9 | 0.4 |

Inter-assay precision

| | Mean | SD | CV |
|----------|--------|------|------|
| n = 20 | U/l | U/l | % |
| sample 1 | 150.03 | 2.36 | 1.57 |
| sample 2 | 484.69 | 6.57 | 1.36 |

5. Methods Comparison

Comparison was done between reference Creatine Kinase Reagent and CliniQuant - FSR Creatine Kinase Reagent (test).

N = 23 y = 0.901x + 38.43

r² = 0.849

LIMITATIONS

Samples with values above 1800 U/l should be diluted with 0.9% saline, re-run and results multiplied by dilution factor.

WASTE DISPOSAL

This product is made to be used in professional laboratories. Please consult local regulations for correct waste disposal.

REFERENCES

- Burtis, C.A., Ashwood, E.R., editors. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia, W.B. Saunders Company, 1994, p. 797 - 809.
- Data on file: Meril Diagnostics.

IFU/CKNF SR01/00

05-12-2018

Symbols used on Meril Diagnostics labels:

| | | | |
|--|--|---|-----------------------------------|
|  REF | Catalogue No. |  | Attention See Instruction for Use |
|  LOT | Batch No. |  | In vitro Diagnostics |
|  | Expiry Date |  | Consult Instruction for Use |
|  | Manufacturer |  | Storage Temperature |
|  | Keep Dry |  | Keep Away from Sunlight |
|  | Manufacturing Date |  | Do not use if package is damaged |
|  | Authorized European Representative in the European Community | | |