

Meril

Diagnostics

ISO 13485 & GMP Certified

Better POWER OF 2

QUADRO™ HIV 1-2 Ab

4 Dot Rapid Test Device with Flow Through Technology for Detection of Anti-HIV Antibody

Evaluated
by
NARI

SENSITIVITY 100%
NARI
SPECIFICITY 99.5%



★
Power of 2
★

- ★ Synthetic Peptides coated on HIV 1 and HIV 2 individual Test Region
- ★ Uniquely formulated Recombinant Proteins Coated on HIV 1 + 2 Test Region

★
Be Extra Sure 

- ★ Synthetic Peptides detects antibodies against the serotype specific linear epitopes
- ★ Recombinant Peptides detects antibodies against the conformational epitopes
- ★ Differential detection of HIV 1 and HIV 2
- ★ Detects antibodies against all group of HIV (M, N and O) Including subtype 'C'
- ★ Sensitivity : 100% and Specificity : 99.5%
- ★ No Hook effect
- ★ Result in **3 minutes**
- ★ User Friendly Protocol



Meril Diagnostics Pvt. Ltd.

D1 - D3, Second Floor, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi - 396191, India

601, Midas, Sahar Plaza Complex, J.B Nagar, Andheri East, Mumbai - 400059

Tel: +91 22 30732400 | Fax: +91 22 39350777 | Customer Care: 1-800-266-3745 (Toll Free) | Email: diagnostics@merillife.com


QUADRO™ HIV 1-2 Ab

4 Dot Rapid Test Device with Flow Through Technology for Detection of Anti-HIV Antibodies

TESTING PROCEDURE



RESULT INTERPRETATION



ORDERING INFORMATION

Product Code	Product Name	Pack Size
HD4RPD-01	Quadro HIV 1 - 2 Ab	10 Tests
HD4RPD-03	Quadro HIV 1 - 2 Ab	50 Tests



Meril Diagnostics Pvt. Ltd.

D1 - D3, Second Floor, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi - 396191, India

601, Midas, Sahar Plaza Complex, J.B Nagar, Andheri East, Mumbai - 400059

Tel: +91 22 30732400 | Fax: +91 22 39350777 | Customer Care: 1-800-266-3745 (Toll Free) | Email: diagnostics@merillife.com